


SUPRAMID

BLACK



| | |
|------------------------|--|
| Type | Twisted and coated multifilament Sizes 5/0 and 6/0 monofilament |
| Composition | Polyamide 6.6 strands enclosed in a polyamide 6 sheet |
| Colour | Black |
| Tissue reaction | Minimal |
| Absorption | Non absorbable, gradually encapsulated by connective tissue |
| Presentation | Box with 12 sutures |

*For quality
&
safety*

smi

www.sutures.be

SUPRAMID

BLACK

CHARACTERISTICS

- ✓ High tensile strength, non absorbable
- ✓ Synthetic
- ✓ Twisted structure

- ✓ Multifilament in a sheath

- ✓ Pliability and softness

- ✓ Hermetically sealed packing

INDICATIONS General soft tissue closing
 and/or ligation
especially in Skin closure

BENEFITS

- Excellent and permanent tissue support

- Minimal tissue reaction

- Excellent handling properties
High flexibility
High tensile strength
Excellent knot security

- Soft passage through the tissues without sawing,
tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action

- Good handling properties
Easy and secure knotting

- Guaranteed seal and product sterility