

# CERVIX SET



## Product

This special suture set consists of a white polyester tape equipped at each end with blunt round bodied needles.

## History

The Shirodkar technique was first described in India in 1955. The Shirodkar stitch is a purse-string stitch put into the cervix to keep it tightly closed during pregnancy.

## Process


A ligature is placed around the cervix in a collar-like fashion, drawn tight, and either sutured together or tied closed. A strong tape is then woven carefully with a swage blunt needle in and out of the mucosa.

## Advantages

- ✓ Will not cut or damage the wall of the cervix, when placed properly
- ✓ Prevents a miscarriage or premature birth when the cervix has been weakened
- ✓ Reduces the danger of premature labor and the associated triggering of a late abortion or a premature birth

## Technical description

Braided non absorbable polyester tape, 4 mm  
12 pc in a box

2 x blunt needle	Product	Order no
	2 x HRB 50 mm	4 mm polyester tape 40 cm <b>ZCE1</b>
		4 mm polyester tape 50 cm <b>ZCE4</b>
	2 x HRB 65 mm	4 mm polyester tape 50 cm <b>ZCE3</b>



*For quality  
&  
safety*



[www.sutures.be](http://www.sutures.be)