

# EPISIO SET



## Product

This special set contains all sutures for a episiotomy procedures in one single foil pack.

## Benefits

- Easy and quick identification of material
- All in one needle suture package

## Features

- Supple passage through tissue
- Fast absorption
- High tensile strength

## Technical description

Surgicryl Rapid: synthetic absorbable Polyglycolic Acid - braided and coated multifilament  
12 pc in a box

### Reverse cutting

Round bodied taper cutting point



Order no

#### 3/8 Circle



DS 26 mm  
HRT 40 mm

1 x USP 2/0 – DS 26 mm – 75 cm  
2 x USP 0 – HRT 40 mm – 90 cm

S-1430-1435

#### 1/2 Circle



1 x USP 2/0 – DS 26 mm – 75 cm  
1 x USP 0 – HRT 40 mm – 90 cm

S-1430-1435-1



For quality  
&  
safety

**smi**  
www.sutures.be