


SMI BIOPSY PUNCH



Product

SMI biopsy punches make it possible to perform skin biopsies in various locations with great precision and minimal tissue trauma

Outstanding quality

- ✓ It's lightweight with a hollow circular, constructed from 304 stainless steel with a seamless razor-sharp cutting edge
- ✓ Cleanly defined length and diameter of the skin sample

Safe collection

- ✓ Ribbed handle provides a secure grip and better control
- ✓ Safety capped for your protection and to maintain sharpness
- ✓ Disposable sterile punch – gamma radiation– single use

Easy to use

- ✓ Size embossed on each punch
- ✓ Single handed ejection cap
- ✓ Sterile & ready to use

Product	Unit	Order no
Biopsy punch Ø 2 mm	25 pc.	ZBP2
Biopsy punch Ø 3 mm	25 pc.	ZBP3
Biopsy punch Ø 4 mm	25 pc.	ZBP4
Biopsy punch Ø 5 mm	25 pc.	ZBP5
Biopsy punch Ø 6 mm	25 pc.	ZBP6
Biopsy punch Ø 8 mm	25 pc.	ZBP8
Biopsy punch Ø 10 mm	20 pc.	ZBP10
Biopsy punch Ø 12 mm	20 pc.	ZBP12



*For quality
&
safety*

smi
www.sutures.be